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# The Implant System for Spine Fixation- Pedicular. BM was projected to use in the Spine by via posterior. It is a Group of metallic implant dystem for given the standard reduceds: Downess projected to use in the cynic by Map posterior. In its a droup of metanic implantable and modulates components, compatible amongst themselves and anatomical to the bone, used for: Correction, Stabilization, and Fixation permanently or temporarily the segments thoracic, lumbar and sacral of the spine. Having Pedicular Poliaxial screw (articulated), Pedicular Screw, Pedicular Spondyloilsthesis Screw, Distraction Stem, Transverse Stem, Connector, Poliaxial screw (articulated), Pedicular Screw, Pedicular spondyloismesis Screw, Distraction Stem, Transverse Stem, Connector, Screw Locks Hook, Hook Normal Distractor, Hook Left Distractor Straight, Hook Lateral Left Distractor, Hook Distractor Soping Right, Hook Distractor, Hook Distractor, Hook Distractor Pedicular and Contra. Screw. All the implants are supplied separately. The finish of the implants is by mechanical polishing afterwards this product goes by a residual removing process, that consists of chemical attack of the polished surface of the metal that removes the residual sludge's of the production process, such

other production materials afterwards this product goes by an ionization process, in other words the superficial layer of the products manufactured in titanium.

### DPODLICTION MATERIAL

PRODUCT DESCRIPTION

The products are made through the titanium alloy 6AI 4V ELI NBR ISO 5832-3 or ASTM F-136. It is a biocompatible alloy and it possesses resistance to the corrosion. The resistance to the fatigue and high resistance to the corrosion is intrinsic mechanical properties of 6AI 4V ELI

# PEDICULAR FIXATION SYSTEM

### ANCILL ARY COMPONENTS

for Spine Fixation-Pedicular BM doesn't possess ancillary compon

### ACCESSORY

ACCESSORY
The Implant System for Spine Fixation- Pedicular BM doesn't possess any accessory with the purpose of integrating the medic product, granting to this product a function or complementally technical characteristic.

MATERIAL OF SUPPORT The support materials that accompany the implants of Implant System for Spine Fixation-Pedicular BM are:
-Bull of the product. Use Instructions. SIFPC (revision 0);

- Surgical technique (optional). (Revision 0); - Box of surgical instrumental (optional). COD.:7000

## FORM OF PRESENTATION

The Implant System for Spine Fixation- Pedicular BM, is wrapped in plastic envelope, identified with no sterile label, external labeling, containing 1 implant unit, 1 instruction of use and 1 internal label to be labeled in the patient's medical file. The hospital institution is responsible in assuring the identification and the traceability of the medical product, through the data and the necessary information, contained in the labeling (label) of the medical product. The labeling contains the following information

Manufacturer Data | Description of packing content | Measure | Raw material | | Init | Reference | No Sterile | Sterilization Method | maintautuel Paia pessipiini ni pacang orinent in measure i naw material print (Peterberge Production Date) Validity [Medical Article of unique use | Instructions, Warnings and Precautions: See instructions o use | Technical responsible | Health Ministery Registration no. | Brazzilian Industry | Product of Unique Use Do not reuse

INDICATION OF USE
The Implant System for Spine Fixation- Pedicular BM, is indicated in the treatment of Lumbar and Thoracic Fractures, by Via
Posterior , disk hemia, pseudarthrosis , Stenosis, Spondylolisthesis, sequels of fractures, resection of Tumors, cifosis, Idiopathic
and congenital Scoliosis and neuro-muscular deformities, Lordosis, Degeneration of Intervertebral Disc and surgical revisions.
NOTE: The implants here indicated they were specifically projected for use above described. Any other type of use is contra

CONTRA-INDICATIONS
It is contra inclinated the use of this implant in the following situations:
-sensibility to the metal; / -fever or sign of local inflammation; / -osteoporosis; / -obesity (load in the excessive spinal system produced by obese patient they can compromise the fixation of the implant for column); / -inadequate tissue covering in the area to be operated; / -pregnancy; / - bone compromising for diseases and/or infections; / -elevation of the tax of sedimentation unexplainable through other diseases, elevation of the tax of WBC\*, or outstanding suspects in differential counting; / -visible intention; / -conscious tracerbuteties, elevation of the tax of WBC\*, or outstanding suspects in differential counting; / -visible intention; / -conscious tracerbuteties, elevation of the tax of WBC\*, or outstanding suspects in differential counting; / -visible intention; / -conscious tracerbuteties, elevation of the tax of WBC\*, or outstanding suspects in differential counting; / -visible intention. nfection; /-precocious deambulation /-other situations, that can disable the success of spinal implant surgery. WBC: White blood counting

### FORESEEN PERFORMANCE IN SAFETY'S ESSENTIAL REQUIREMENTS AND EFFECTIVENESS OF THE PRODUCT.

Modeling Limits or Mechanical Conformation of Metallic Implantable Components

For the success of the surgery it is fundamental the correct selection of the implants. It is of extreme importance the correct handling of the implants. When accomplishing the modeling in the implants, should be used appropriate equipment. The Distraction Stem and the Transverse Stem were projected to adapt appropriately to the spine segment. During the modeling it should be avoided scratches, deep marks in the implants because those defects can provoke internal tensions in the implant and could turn to a ossible fracture focus in the implant. In case of necessity this components can be modeled or conforme

possible fracture focus in the implant. In case of necessity this components can be modeled or conformed mechanically respecting as modeling limit the own bone anatomy, avoiding folds in acute angles, inverse folds and repetitive modeling. The remodeling or a new mechanical conformation for his original form or inverse fold is prohibited, because the implant will suffer mechanical tensions that will cause the break of the implant. The assembly of the system should be careful allowing the same to be fit in perfectly in the parts that contain threads, mainly in the instant of placement of the obstacle screw in the screw he dicular. This needs to be put carefully for the squeeze to be effective. Damages provoked in the thread during the surgery or if the nut enters pie in the screw the appropriate squeeze it will be committed. The items below describe results of rehearsof load restrictions, limits of mechanical resignation and torque in the screw and accomplished with the Implant System for Spine Fixation-Pedicular BM.

### Restrictions of Load (Summarize of Rehearsals)

strictions of Load (Summarize or Renearsals)
mpression rehearsals: Tests demonstrated that the Implant System for Spine Fixation- Pedicular. BM supporte
1407,37 N: Load Comp. Esc. 1296,00 N; rigidity Compressive. 270,97 N/mm; Deviation of Desl. 6,65mm; 4,94mm
The methods for obtaining these limits are established in the norm ASTM F1717-04 ((Standard Test Method for Spinal Implant

Constructs in a Vertebrectomy Model). The rehearsal in the Implant System for Spine Fixation- Pedicular BM was accomplished as specified in the norm ASTM F1717-04 Standard Test Method for Spinal Implant Constructs in a Vertebrectomy Model). Bar sliding: Sliding Rehearsal of bar/screw demonstrated that the System of Implants for Fixation of Spine Pedicular BM supported medium load of 236,71Kgf. The Sliding rehearsals were accomplished using a constant speed of 5mm/min.

Traction: In the traction rehearsal the raw material used for production of the System of Implants for Fixation of Spine Pedicular BM resented limit of drainage of 918 MPa, limit of resistance of 963 Mpa and 17% of prolongation assisting the normASTM F136-02a

### INFORMATION OF USE

or the correct use of the implantable, metallic components of the Implant System for Spine Fixation Pedicular-RM, is necessary For the correct use of the implantable metallic components of the Implant System for Spine Fixation Pedicular- BM, is necessary:

1. Deep anatomical knowledge of each vertebral segment to suffer surgery and knowledge of the shology to be treated: 2. The surgeon's domain of the surgical procedure with osteosinthesis and arthrodesis notions and because the system has as purpose to correct, to stabilize and to latent he segments thoracic, lumbar and sacral of the spine. Finally, the use of the Implant System for Spine Fixation Pedicular- BM should be made by skilled surgeons for this procedure. 3. Analysis of the patient's bone structure verifying if supports the dimensions of the implants of Implant System for Spine Fixation Pedicular- BM. 4.1otal care is necessary in the access, in the perforation of the pedicles and entrance of the screws. We recommended for a better safety the use of intensifier of images, because during the perforation of the pedicles and entrance of the screws all care is necessary. 5. The hospital should fill out the patient's records scoring the code of the product, number of the lot of the implant. Biomecanica sends 1 adhesive label to be glued in the patient's records scoring the code of the product, number of the lot of the implant. Biomecanica sends 1 adhesive label to be glued in the patient's records coming the cord all the data of the implant mentioned above. 6. It is important the use of the surgical instrumental was specifically projected for use with this implant system. 7. The implants of the Implant System for Spine Fixation Pedicular- BM are of unique use. The resultization is EXPRESSLY PROHIBITED.

A) Sensibility reactions to the metal in patients were rarely informed. Implantation of strange material in tissues results in histological reactions that involve several sizes of macrophages and fibroblasts. The chemical importance of this effect is uncertain, as well as similar changes can happen as a precursor, or during the healing process; b) Long bone coolidation that can take to the break of the implant; c) Mechanical loosening, sprain, dismemberment, breaks of the components - they can be the result of defective fixation or infection hides; d) infection superficial and/or deep; e) Shortening of the segment and of the spine in function of bone reabsorption; pain, discomfort or abnormal sensations due to placement of the implantable component metallic or due to migration of the implantable mechanical components thrombosis; f) Bone fractures during the placement of screws and hooks; g) Pseudo-arthritis; h) Fracture; i) Paralysis; j) bruises and hemorrhage in the blood vessels; l) lung clot; m) neural problems or neurological caused by the surgical trauma as: gastro-intestinal disturbances; n) - Precocious or late fails, of the implantable mponents ; o) scars; p) superficial or muscular sensibility in patients with inadequate of tissue covering in the operated area; q) difficulty in accomplish physical activities; r) allergy; s) Illness and Death; in all of the

surgical procedures exists an incidence of illness and death. The patient should be communicated by the surgeon before suffering

### PRECAUTIONS

PRECAUTIONS

The surgeon should be aware of the possibility of the patient's physical development, could have need of revision of the implant;
The surgeon should be familiarized and to stop enough knowledge of the surgical technique of spinal fixation and his limitation, including the pre and postoperative, adopted surgical technique, precautions and potential risks;
The surgeon should have attention in the preparation of the pedicile for the introduction of the pediculars screws and modeling of the bars, because the stability and the success of the fixation can be committed;
It is of extreme importance the correct handling of the implants. When accomplishing the modeling in the implants, appropriate equipment should be used. Scratches, deep marks and inverse fold should be avoided in the implants during the modeling, because those defects can provoke internal tensions in the implant could turn a possible fracture focus in the implant;
To guarantee an appropriate implantation, the instrumental of column should be used manufactured by the Biomecanica, because these were projected and specifically manufactured for the use of those implants. Use instrumental of different manufacturers it can commit the surgery;

ommit the surgery; The implants cannot be bent in an accentuated way or inverted;

-Not to use the product if it is damaged; -PRODUCT NO STERILE. STERILIZE BEFORE USE; -PRODUCT OF UNIQUE USE. IT IS EXPRESSLY PROHIBITED THE RE-USE OF THIS PRODUCT;

Special conditions of storage: Conserve at airy place, dry and to the shelter of the light and far away from the action of bad · Do not use the product in case the packing is damaged and/or period of production expired: Observation: Components of different manufacturers' implants should not be used: therefore we recommend that the products

### THE PATIENT SHOULD BE INFORMED AS:

To this instruction in the items: Indications, Contra-Indications, Information of Use - Possible Adverse Effects, Precautions and Warnings,

nave the same origin

The patient should be informed as the importance of the post surgical attendance. The lack of attendance impedes the detection of

Phone: 0034 91 517 2492 Fax: 0034 91 842 9180

MANUFACTURER'S AUTHORIZED EU REPRESENTATIVE

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### WADNINGS

roduct of Unique Use. Destroy after explanted. Not to reuse the product. An orthopedic implant can only be used in a single patient, a single time. Although it can seem not to be damaged, previous tensions can create imperfections that can reduce the success of the implant. Inappropriate selection of implant can cause uncommon tensions in the implant and they can result in subsequent fracture of the implant.

and they arressum subsequent relative memplant.

Due to his presentation form and their characteristics, this product doesn't present side effect, however, we suggests monitoring the patient's re-establishment through the X-Ray, according to International Norms of Safety, Attention when using magnetic resonance in patients that already has vocume implants manufactured in stainlessel, because it erence in the exam harming the diagnoses.

cause interierce in the examinarming the diagnoses. The use of the Implant System for Spine Fixation Pedicular BM is conceived for correction, stabilization and fixation of the segments thoracic, lumbar and sacral of the spine and not to substitute the normal structures of the human Obesity has effect in the load in the spinal system could commit the fixation of the implant for column. Excessive physical

activity and traumatisms called directly in the segments thoracic, lumbar and sacral can commit the success of the implant prematurely. We recommended that the activities are dosed avoiding the excessive "stress" in the operated area. The n for Spine Fixation Pedicular, RM cannot suffer loads or levels of similar activities supply

The mixture of different biomedical materials no compatible and/or similar it can elevate the corrosion risk taking int account also that in spite of the implant to be manufactured by material biocompatible when implants in the patient the same is subject to constant environmental changes promoting the concentration of salts, acids and alkaliotis that it ca promote the corrosion. The clinical results and durability of the implant depend on the existence of a surgical technique. needs. He/she dates from production, period of validity and lot of the product: VIDE LABEL.

### INSTRUCTION OF USE

INSTRUCTION OF USE
The Surgical Techniques depend on the technique adopted by the surgeon, fitting to this the choice of the method, type and dimension of the products to be inserted.
The mixture of different biomedical materials no compatible and/or similar it can elevate the corrosion risk. The clinical results and durability of the implant depend on the existence of a surgical technique needs.

The placement of the Implant System for Spine Fixation- BM is accomplished with aid of the Instrumental Kit of Column Biomecanica

### Preparation for screws placemen

I Use the Perforator 8509 for rip off the point of entrance of the pedicle, exposing the spongy structure. The point of entrance of the pedicle is intersected by the vertical line which connects the lateral borders of the bone crest extension of intrarticular segments, and the horizontal line which divides the transverse process. Controlling the apide pressure in intrarticular segments, and the nonzontal nine wind individes the transverse process. 2. Controlling the applied pressure in a series of rotating and firm movements, the meter 8510 is inserted in the pedicle with angulation from 10 to 20°. This meter will follow the timerary of smaller resistance (spongy structures). 3. The feeler guide 8514 is used to verify the presence of bone plan in all of the directions. The tip of the instrument should slide in contrary sense to the cortical wall in all the plans, 4. In case of doubts, X ray should be accomplished with help of the meter pin 8513.

Insert of the screw

5.Use the square reamer tip 8511 to determine the appropriate length of the screw. The extremity of pin depth meter is put in the channel of the pedicle and it is advanced through the channel until the previous wall of the vertebral body, 6. The screws are inserted using the Initial Key Guide 8524. The tip of the key is inserted in the compartment of the pedicular screws are inserted using the Initial Key Guide 8524. The tip of the key is inserted in the compartment of the pedicular screw and tight locking the system. O88: The pedicular screws are auto-reaming, usually doern't request thread preparation. However, the males for this preparation are available in the Instrumental Kit of Column.-Biomecanica. It is important to remind that the pedicular screws should not be removed for the depth adjustment. In consequence, should introduce each aligned screw in relation to the adjacent ones. This way checks that the distraction stem will follow a regular line in the correct form.

Insert of the distraction stem

7. After placement of pedicular screws, the distraction stem is selected. The same can be modeled with the aid of Stem Modeler 8518, when necessary. The stems are supplied in varied lengths and they should be extended 5 mm approximately besides the external borders of the screws bodies of superior pedicle, and more inferior. 8. The Fillister Stem Guide 8502 can be used to stabilize the position of the Distraction Stem as the blockade screw of pedicular screw was not inserted.

# was not inserted. Insert of transverse stem

9. In certain situations the transverse stem should be applied to provide rigidity to the system and to distribute the loads better among the vertebrae pedicles. The size of the transverse stem should be defined so his compression is appropriated, as well as it can promote the appropriate fixation of the whole group. 10. Use the transverse hook coupled to the transverse stem to establish the fixation of the group.

### RISKS OF CONTAMINATION

Risks of biological contamination and transmission of viral diseases exist as HIV and hepatitis because the metallic components of the Implant System for Spine Fixation Pedicular - BM gets in touch with tissue and coproral fluids. Products explanted should be treated as highly pollutants.

### METHODS AND PROCEDURES OF THE MEDICAL PRODUCT DISCARD

The metallic components of Implant System for Spine Fixation Pedicular - BM explanted of patients should be discarded properly by the hospitalar institution. Is of entire hospitalar institution responsibility. the method to discard the implant impeding his/her reverse-use. Biomecanica recommends that the explanted implants is deformed mechanically with aid of hammer, impact press, should be soon afterwards identified with the words "Inappropriate for Use."

This product is supplied no sterile. Before use should be sterilized. We recommended the steam sterilization in h Autoclave (ISO 11134: 1994 Sterilization of health care products Requirements for validation and routine Industrial moist heat sterilization).

Autoclave
It is equipment for sterilization through saturated steam under pressure.
Picture: Sterilization of the Material for the steam process saturated under pressure (Autoclave), In agreement with the temperature exposition according to (Resolution SS-374, de 15-12-95) Time of Exposition

rempe	erature	Time of Exposition				1:
Materi	ial	120° C	130°C	132° C (High Values)	Observation	i
	al of stainless steel er metal thype	30min	15min	4min	Wrapped in paper kraft or field of fabric of double raw cotton.	
	ic instruments inserted ys or metallic box ated	30min	15min	4min	Wrapped in paper kraft or field of fabric of double raw cotton.  Contra indicated the use of boxes semi open	E U
Cutting and sa	g blades scissors ww	30min	15min	4min	Wrapped up individually or conditioned in tray or perforated and wrapped up metallic box in paper kraft or field of fabric of double raw cotton.	1

FOR BETTER EXPLANATION CONSULT THE MANUAL OF INSTRUCTIONS THAT ACCOMPANIES EACH

**Operational instructions** For departure of the equipment, verify: if the interrupter is linked; / if the water register is open; / if the discharge register

closed.

1. Open the door of the equipment; 2. Accommodate the material to be sterilized appropriately; 3. Close the door of the equipment; 4. Select the cycle wanted in agreement with the material to be sterilized; 5. Tie the general key; 6. Cycle will elapse automatically, in the sequence; 7. When lighting the "final lamp of cycle", open the door partially for ten (10) minutes, approximately, for cooling the material.

# FOR BETTER EXPLANATION CONSULT THE MANUAL OF INSTRUCTIONS WHICH ACCOMPANIES EACH

Another sterilization method: besides AUTOCLAVE can be used the following method of sterilization in accordance defined by hospitalar institution:

Sterilization by oxide of ethylene (E.T.O.) - parameters and procedures established in the validation protocol and in FN550:1994 Sterilization of Medical Devices Validation and Routine Control of Ethylene Oxide Sterilizati NOTE

It is of entire care the hospitalar institution responsibility, for the sterilization method, the equipments, the controls, and the

# CLEANING AND DECONTAMINATION RECOMMENDED

C€

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General recommendations
To guarantee the reduction of the infection risks, the implants and surgical instruments must be cleaned after the st and before the sterilization. The products should have decontamination using a bactericidal solution and fungicide o spectrum, avoiding aggressive agents of cleaning, metal brushes, so the quality of the metallic implaits is not han The manual cleaning is more recommended by attacking less the implants and surgical instruments. In the m cleaning, water should be used at room temperature, t use nylon brushes, never of steel or abrasive, because it can the protecting layer of the material.

BIO MECANICA

In case of cleaning machines the implants and surgical instruments should be positioned in drawers, to avoid hit one in another, avoiding damages in the material. It is not recommended to put together materials of different composition luring the wash in cleaning machines f

The responsible for cleaning the implants and surgical instruments should be attentive as for the type of cleaning product used respecting the suitable concentrations by the manufacturer. The implants and surgical inst be rinsed carefully until every residue comes out dry immediately after the cleaning

# RE-STERII IZATION

ion, all the cares above described should be assisted

### TRACEABILITY

Near the packing of the implantable component, follows one more label containing the data of the used implant. This abel should be glued in the patient's file. The hospital should register that information in their files and to review that patient information. In the file the following information are indispensable: Name of the implant/- Surgery dates / - code of the product/- lot number/- patient's name / - surgeon's name

All the implants mentioned in this instruction of use, count the following information recorded by laser to allow the aceability and the postoperative patient's attendance:
Manufacturer's logo /- size /- lot /- reference

### SPECIAL CARES IN THE STORAGE AND TRANSPORT OF THE MEDICAL PRODUCT

- Do not use the product if he looks to be damaged. / - The implants should be manipulated carefully. Avoiding damages that can harm the quality of the material and the patient's safety /-Conserve at airy place, dry, to the shelter of the light and far away from the action of bad weather. / - When transporting t, shocks, vibrations and high temperatures shou be avoided (above 45°) and inadequate piling up.

The cares to be taken during the transport are described in the collective packing.

### OBSERVATION

This product should only be manipulated by qualified people omponents of different manufacturers' implants should not be used. Therefore we recommended the products have

Every implant is of unique use. Therefore, when withdraw from patient, should be clean and discarded

### COMPLAINT

COMPLANTAI
In case of some complaint of Implant System for Spine Fixation BM these should be reported to Biomecanica authorized distributor. In case the complaint is related to some adverse effect that affects the user's safety, as producing working, harm of the metallic component, serious problems or death related with those components Biomecanical

SENDING THE MATERIAL FOR MANUFACTURER ANALYSIS
In case of sending implants for manufacturer accomplish analysis, this should be higienized in the hospital using a bactericidal solution and fungicide of wide spectrum. Afterwards it should be disinfected or sterilized by oxide of ethylene. They should be dispatch to Biomecanica in integrate packing, identified with the method of cleaning, terilization and data of the product.

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